


## KRAS, BRAF, PIK3CA Array

*Targeted mutation profiling*



# KRAS, BRAF, PIK3CA\* Array


## Rapid profiling of point mutations in the KRAS, BRAF and PIK3CA genes

### Introduction

Colorectal cancer (CRC) is the third most common cancer worldwide.<sup>1</sup> Metastatic disease accounts for 40-50% of newly diagnosed patients and is associated with high morbidity.<sup>1,2</sup> Despite recent therapeutic advances, the prognosis for patients with metastatic CRC (mCRC) remains poor.<sup>3</sup> In recent years monoclonal antibodies (moAbs), like cetuximab and panitumumab which target the epidermal growth factor receptor (EGFR), have proven to be effective in combination with chemotherapy or as single agents for the treatment of mCRC.<sup>4,5</sup> These moAbs

block the signal from EGFR inhibiting downstream signalling including KRAS, BRAF and PIK3CA mediated events (see diagram below). However, when KRAS, BRAF and PIK3CA\* are mutated they are permanently 'turned on' permitting downstream events irrespective of anti-EGFR therapy. The Randox KRAS, BRAF, PIK3CA Array allows the clinician to detect important mutations in the KRAS, BRAF and PIK3CA genes, enabling the appropriate selection of patients for therapy.

### EGFR pathway and its inhibition by EGFR targeted monoclonal antibodies



### Targets detectable by the KRAS, BRAF, PIK3CA\* Array

The KRAS, BRAF, PIK3CA\* Array is designed for the rapid qualitative detection of point mutations within the genes KRAS, BRAF and PIK3CA from fresh/frozen and formalin fixed paraffin embedded (FFPE) tissue DNA (refer to below table for targets detectable).

KRAS						BRAF	PIK3CA*		
codon 12		codon 13	codon 61		codon 146	codon 600	codon 542	codon 545	codon 1047
G12A	G12C	G13D	Q61K	Q61H(1)	A146T	V600E	E542K	E545K	H1047R
G12R	G12S	G13C	Q61L	Q61H(2)	A146P				
G12D	G12V	G13R	Q61R						

\*PIK3CA for research use only

## Why test the KRAS, BRAF, and PIK3CA genes?

Early studies conducted on mainly heavily pre-treated chemotherapy-refractory patients and also chemotherapy-naïve patients with mCRC indicated that only 10-20% of patients clinically benefited from anti-EGFR moAbs.<sup>5,6</sup> Consequently oncogenic activation of EGFR downstream effectors was investigated with respect to clinical outcome to moAb therapy. Analysis confirmed that patients with mCRC carrying activating KRAS gene mutations do not benefit from anti-EGFR moAb therapy.<sup>7,8</sup> KRAS mutations have since emerged as the major negative predictor of

efficacy in patients receiving cetuximab or panitumumab.<sup>3</sup> The occurrence of KRAS mutations however only accounts for approximately 35-45% of nonresponsive patients.<sup>3</sup> The identification of additional genetic determinants of primary resistance to EGFR-targeted therapies in colorectal cancers is therefore important. Recent studies have focused on the molecular analysis of the molecules involved in downstream EGFR signalling with mutations in BRAF<sup>9</sup> and PIK3CA<sup>10</sup> genes being reported to affect patient response to EGFR-targeted moAbs.

## KRAS, BRAF, PIK3CA\* Array Protocol



## Benefits of the KRAS, BRAF, PIK3CA Array

### To the patient

- KRAS, BRAF, PIK3CA\* Array is a rapid simple method for mutation detection
- Determining mutational status informs selection of appropriate therapy

### To the laboratory

- Detection of 1% mutant in a background of wildtype genomic DNA
- Compatible with a broad range of genomic DNA input and type:
  - Formalin fixed paraffin embedded (FFPE) tissue
  - Fresh/frozen tissue
- Single DNA sample required
- Turnaround time of 3 hours
- Single reaction multiplex PCR coupled to a biochip provides greater mutation coverage of the three most important genes (KRAS, BRAF and PIK3CA) implicated in metastatic colorectal cancer therapy response
- Streamlined workflow – protocol and reagents are optimised for the molecular laboratory

# Evidence Investigator

Multiplexing...proven, perfected, evolved

The Evidence Investigator is a semi-automated, benchtop biochip analyser which offers complete patient profiling.

## Save time and costs

Multiplexing reduces time, labour and reagents associated with multiple individual tests

## Increase throughput

For greater laboratory efficiency

## Consolidation

Of immunoassays and molecular diagnostics, improving laboratory efficiency

## Result traceability

Chain of custody features and bar coded reagents

## No hidden costs

Package includes imaging module, PC and imaging software, thermoshaker, biochip carrier handling tray and barcode scanner

## Ease of operation

Straightforward testing procedure, ready-to-use biochips and minimal sample handling

## Extensive QC

Internal quality controls ensure all key assay steps have been performed correctly i.e. amplification

## Retrospective reporting

Enabling additional analysis of previously captured sample data

## Ordering Details

Description	Size	Cat. No.
KRAS, BRAF, PIK3CA Array	54 Biochips	EV3799A & EV3799B
Evidence Investigator Analyser		EV3602

\*Note: Extraction reagents are not included  
**PIK3CA\* for research use only**

## References

1. Ferlay, J., Autier, P., Boniol, M., et al. (2007) Estimates of the cancer incidence and mortality in Europe in 2006. *Ann Oncol* 18, 581-592.
2. Centers for Disease Control and Prevention: United States Cancer Statistics: US Cancer Statistics Working Group. <http://www.cdc.gov/uscs>.
3. Bardelli, A. & Sienna, S. (2010) Molecular Mechanisms of resistance to cetuximab and panitumumab in colorectal cancer. *Journal of Clinical Oncology* 28(7) 1254-1261.
4. Salt, L.B., Meropol, N.J., Loehrer, P.J. Sr. et al. (2004) Phase II trial of cetuximab in patients with refractory colorectal cancer that expressed the epidermal growth factor receptor. *Journal of Clinical Oncology* 22, 1201-1208.
5. Cunningham, D., Humblet, Y., Siena, S., et al. (2008) Cetuximab monotherapy and cetuximab plus irinotecan in irinotecan-refractory metastatic colorectal cancer. *New England Journal of Medicine* 358, 1160-1174.
6. Chung, K.Y., Shia, J., Kemeny, N.E. et al. (2005) Cetuximab shows activity in colorectal cancer patients with tumours that do not express the epidermal growth factor receptor by immunohistochemistry. *Journal of Clinical Oncology* 23, 1803-1810.
7. Lièvre, A., Bachet, J.B., Le Corre, D. et al. (2006) KRAS mutation status is predictive of response to cetuximab therapy in colorectal cancer. *Cancer Research* 66, 3992-3995.
8. Benvenuti, S., Sartore-Bianchi, A., Di Nicolantonio, F. et al. (2007) Oncogenic activation of the RAS/RAF signalling pathway impairs the response of metastatic colorectal cancers to anti-epidermal growth factor receptor antibody therapies. *Cancer Research* 67, 2643-2648.
9. Di Nicolantonio, F., Martini, M., Molinari, F. et al. (2008) BRAF wildtype is required for response to panitumumab or cetuximab in metastatic colorectal cancer. *Journal of Clinical Oncology* 26, 5705-5712.
10. De Roock, W., Claes, B., Bernasconi, D. et al. (2010) Effects of KRAS, BRAF, NRAS and PIK3CA mutations on the efficacy of cetuximab plus chemotherapy in chemotherapy-refractory metastatic colorectal cancer: a retrospective consortium analysis. *Lancet Oncology* 11, 753-762.



**RANDOX**  
MOLECULAR

Randex Laboratories Limited, 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom  
T +44 (0) 28 9442 2413 F +44 (0) 28 9445 2912 E [marketing@randox.com](mailto:marketing@randox.com) I [www.randox.com](http://www.randox.com)

Information correct at time of print.  
Randox Laboratories Limited is a subsidiary of Randox Holdings Limited a company registered within Northern Ireland with company number NI. 614690. VAT Registered Number: GB 151 6827 08. Product availability may vary from country to country. Please contact your local Randox representative for information. Products may be for Research Use Only and not for use in diagnostic procedures in the USA.



LT345 NOV14