RANDOX

CYSTATÍN C

---

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## REFERENCIE

1. Bello, AK, et al. Global Kidney Health Atlas: A report by the Internal Society of Nephrology on the current state of organization and structures for kidney care across the globe. Brussels : Internal Society of Nephrology, 2017.
2. Bikbov, Boris. Chronic kidney disease: impact on the global burden of mortality and morbidity. *The Lancet*. [Online] 2015. <http://www.thelancet.com/campaigns/kidney/updates/chronic-kidney-disease-impact-on-global-burden-of-mortality-and-morbidity>.
3. National Kidney Foundation. Global Facts: About Kidney Disease. National Kidney Foundation. [Online] National Kidney Foundation, 2015. [https://www.kidney.org/kidneydisease/global-facts-about-kidney-disease#\\_ENREF\\_1](https://www.kidney.org/kidneydisease/global-facts-about-kidney-disease#_ENREF_1).
4. Neuen, Brendon Lange, et al. Chronic kidney disease and the global NCDs agenda. s.l. : BMJ Global Health, 2017.
5. National Institute for Health and Care Excellence. Chronic kidney disease in adults: assessment and management: 1 Recommendations. National Institute for Health and Care Excellence. [Online] January 2015. <https://www.nice.org.uk/guidance/cg182/chapter/1-recommendations#classification-of-chronic-kidney-disease-2>.
6. The Renal Association. CKD stages. The Renal Association. [Online] no date. <https://renal.org/information-resources/the-uk-eckd-guide/ckd-stages/>.
7. Lascano, Martin E and Poggio, Emilio D. Kidney Function Assessment by Creatinine-Based Estimation Equations. Cleveland Clinic. [Online] August 2010. [Cited: 16 May 2018.] <http://www.clevelandclinicmeded.com/medicalpubs/diseasemanagement/nephrology/kidney-function/>.
8. Swedko, Peter J, et al. Serum Creatinine Is an Inadequate Screening Test for Renal Failure in Elderly Patients. Research Gate. [Online] February 2003. [Cited: 6 May 2018.] [https://www.researchgate.net/publication/8243393\\_Serum\\_Creatinine\\_Is\\_an\\_Inadequate\\_Screening\\_Test\\_for\\_Renal\\_Failure\\_in\\_Elderly\\_Patients](https://www.researchgate.net/publication/8243393_Serum_Creatinine_Is_an_Inadequate_Screening_Test_for_Renal_Failure_in_Elderly_Patients).
9. Fassett, Robert G, et al. Biomarkers in chronic kidney disease: a review. Science Direct. [Online] 2 October 2011. [Cited: 16 May 2018.] <https://www.sciencedirect.com/science/article/pii/S0085253815551399>.
10. Mishra, Umashankar. New technique developed to detect chronic kidney disease. Business Line. [Online] 07 May 2018. [Cited: 17 May 2018.] <https://www.thehindubusinessline.com/news/science/new-technique-to-detect-chronic-kidney-disease/article23803316.ece>.
11. Chew, Janice SC, et al. Cystatin C-A Paradigm of Evidence Based Laboratory Medicine. NCBI. [Online] 29 May 2008. <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2533150/>.
12. Zhang, M, et al. Clinical evaluation of serum cystatin C and creatinint in patients with chronic kidney disease: a meta-analysis. NCBI. [Online] 04 August 2013. [Cited: 19 April 2018.] <https://www.ncbi.nlm.nih.gov/pubmed?term=23760917>.
13. Nephrol, Indian J. Serum cystain C as a marker of renal function in detection of early acute kidney injury. NCBI. [Online] 23 May 2013. [Cited: 19 April 2018.] <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3692142/>.
14. National Institute for Health and Care Excellence. Chronic kidney disease in adults: assessment and management: 2 Implementation: getting started. NICE. [Online] January 2015. [Cited: 19 April 2018.] <https://www.nice.org.uk/guidance/cg182/chapter/implementation-getting-started>.
15. Grubb, Anders. Cystatin C is Indispensable for Evaluation of Kidney Disease. NCBI. [Online] 28 December 2017. [Cited: 19 April 2018.] <https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5746836/>.
16. National Kidney Foundation. Cystatin C: What is its role in estimating GFR? National Kidney Foundation. [Online] 23 April 2008. [Cited: 24 April 2018.] [https://www.kidney.org/sites/default/files/02-10-0204\\_GA\\_CystatinC.pdf](https://www.kidney.org/sites/default/files/02-10-0204_GA_CystatinC.pdf).
17. Woo, Kwang-Sook, et al. Clinical Usefulness of Serum Cystatin C as a Marker of Renal Function. NCBI. [Online] 20 Auguts 2014. [Cited: 24 April 2018.]
18. NHS. Chronic kidney disease. HNS. [Online] 15 August 2016. [Cited: 14 April 2018.] <https://www.nhs.uk/conditions/kidney-disease/treatment/#lifestyle-changes>.
19. Herget-Rosenthal, S, et al. Early Detection of Acute Renal Failure by Serum Cystatin C: A New Opportunity for a Hepatologist. Wiley Online Library. [Online] June 2005. [Cited: 19 April 2018.] <https://aasldpubs.onlinelibrary.wiley.com/doi/pdf/10.1002/lt.20422>.

## POZNÁMKY:

A series of horizontal dotted lines for writing notes.

**RANDOX**  
REAGENTS

Copyright © 2019 Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151682708. Product availability may vary from country to country. Some products may be for Research Use Only. For more information on product application and availability, please contact your local Randox Representative.