



ACCURATELY DIAGNOSE VETERINARY PATIENTS WITH ON-SITE TESTING



Introducing the RX misano, a compact, semi-automated chemistry analyser ideal for veterinary testing on a variety of species.

Register your interest at:

randox.com/rx-misano-veterinary



» COMPREHENSIVE, CUSTOMISABLE VETERINARY TEST MENU EXPANDS DIAGNOSTIC CAPABILITIES



Sample volume of 2-20µl makes it ideal for the diagnosis of small and large animals



USB port for importing test menus and exporting results in CSV format allows for easy storage of patient data



Touchscreen technology and icon based software enable ease of use and smooth navigation



Custom build test menus and create favourites to enable personalised panels

VETERINARY TEST MENU

Albumin
Aldolase
Alkaline Phosphatase
ALT (GPT)
Ammonia
Amylase
AST (GOT)
Bile Acids
Bilirubin (Direct)
Bilirubin (Total)
Calcium
Chloride
Cholesterol
Cholinesterase
CK-NAC
CO₂ Total
Copper
Creatinine
D-3 Hydroxybutyrate (Ranbut)

Gamma-GT
GLDH
Glucose Oxidase
Glucose Hexokinase
Glutathione Peroxidase (Ransel)
Glycerol
Iron (UIBC)
Lactate
Lipase
Magnesium
NEFA (Non-esterified Fatty Acids)
Phosphorus (Inorganic)
Potassium
Sodium
Superoxide Dismutase (Ransod)
Total Protein
Triglycerides
Urea
Uric Acid

Urinary Protein
Zinc

Other Tests Also Available:

Acid Phosphatase
Apolipoprotein A-I
Apolipoprotein B
ASO
Beta-2 Microglobulin
CK-MB
Cystatin C
Direct HDL-Cholesterol
Direct LDL Cholesterol
Glutathione Reductase
G-6-PDH
LDH
Microalbumin
Total Antioxidant Status

RANDOX

A global market leader in the in vitro diagnostics industry for over 35 years

+44 (0) 28 9442 2413 • therxseries@randox.com • randox.com/rx-misano-veterinary

Product availability may vary from country to country. Please contact your local Randox representative for information.

Products may be for Research Use Only and not for use in diagnostic procedures in the USA.