

1. Information of the substance/mixture and of the company/undertaking	
1.1 Commercial product name	Lithium 2 (RX Suzuka)
Cat. No.:	LM8053
1.2 Company:	RANDOX Laboratories Ltd., 55 Diamond Road, Crumlin, Co. Antrim, United Kingdom, BT29 4QY Tel: +44 (0) 28 9442 2413; Fax: +44 (0) 28 9445 2912 Email: <a href="mailto:msds@randox.com">msds@randox.com</a> Website: <a href="http://www.randox.com">www.randox.com</a>
1.3 In Emergencies:	Tel: +44 (0) 28 9442 2413

2. Hazards identification	
	Buffer - Contains Potassium Hydroxide (6.6%)
Classification:	Buffer - C: Corrosive

3. Composition / Information on Ingredients	
Warnings:	For <i>in vitro</i> diagnostic use only. Do not pipette by mouth. Handle laboratory reagents in accordance with Good Laboratory Practice. May cause irritation to skin and eyes, may be irritating to mucous membranes and upper respiratory tract. Avoid contact with skin, inhalation and ingestion.
Buffer Risk phrases:	R35; Causes severe burns
Safety phrases:	S24/25; Avoid contact with eyes and skin S26; In case of contact with eyes, rinse immediately with plenty of water and seek medical advice S36/37/39; Wear suitable protective clothing, gloves and eye / face protection S45; In case of accident or you feel unwell seek medical advice immediately (show label where possible).

4. First aid measures	
Emergency First Aid:	Call an internal First Aider where available or else contact a physician immediately. Direction for First Aider or Physician:
Inhalation:	If inhaled, remove to fresh air. If breathing becomes difficult, call a physician.
Ingestion:	If swallowed wash out mouth with water provided person is conscious, call a physician.
Skin Contact:	In case of skin contact immediately wash skin with soap and copious quantities of water.
Eye Contact:	Flush with copious amounts of water for at least 15 minutes. Assure adequate flushing by separating eyelids with fingers. Call a physician.

5. Fire-fighting measures	
Flash point:	Not determined
Flammable limits:	Not determined
Extinguishing media:	No restriction
Special Fire Fighting Procedures:	Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.
Unusual Fire and Explosion Hazards:	In a fire or if heated, a pressure increase will occur and the container may burst.
Thermal decomposition:	Not determined

6. Accidental release measures	
After spillage:	<p>Use appropriate spill absorbent kit as instructed by the manufacturer. Alternatively, mop up with an absorbent material and hold for waste disposal.</p> <p>Appropriate Personal Protective Equipment should be worn, e.g. laboratory coat, gloves, safety glasses and mask.</p> <p>Ventilate area and wash spill site after material pick up is complete.</p>
Absorbent material:	No restriction

7. Handling and storage	
7.1 Handling:	Wash thoroughly after handling. Do not use if skin is cut or scratched. Avoid contact with eyes, skin and clothing.
7.2 Storage	Store at temperatures and conditions as indicated on the product label
Do not store together:	N/A

8. Exposure controls/personal protection	
MEL/OES	Not determined
Respiratory protection:	None
Eye protection:	Protective glasses
Hand protection:	Standard laboratory rubber or latex gloves.
Skin protection:	A laboratory coat or apron is recommended

<b>9. Physical and chemical properties</b>	
Physical state:	Liquid
Colour:	Colourless
Odour:	Odourless
pH value:	Buffer 13.5
Boiling point	Not determined
Melting point:	Not determined
Oxidising properties:	Not determined
Ignition temperature	Not determined
Explosion limits:	Not determined
Density:	Not determined
Solubility in water:	Soluble

<b>10. Stability and reactivity</b>	
Stability:	Stable
Conditions to avoid:	Not determined
Substances to avoid:	Not determined
Hazardous decomposition products:	Under normal conditions of storage and use, hazardous decomposition products should not be produced.

<b>11. Toxicological information:</b>	
LD50 oral	Rat: 273 mg/kg
LD50 intraperitoneal	Not determined
LD50 skin	Not determined
LC50 inhalation	Not determined
Toxicologic information:	Not thoroughly investigated

<b>12. Ecological information</b>	
Water hazard class:	Data not yet available

<b>13. Disposal considerations</b>	<p>Each disposal facility must determine proper disposal methods of the substance or preparation and any contaminated packaging to comply with Local and National Environmental Regulations</p> <p>Refer to Section 6.</p>
------------------------------------	--

14. Transport Information			
GGVE/GGVS:	Data not yet available	RID/ADR	Data not yet available
ICAO/IATA-DGR:	Data not yet available	UN-Nr.:	Data not yet available
IMDG-CODE:	Data not yet available	MFAG	Data not yet available
		EMS:	Data not yet available

15. Regulatory information	
	Classified according to CHIP 3.
Hazard symbol:	Buffer - C; Corrosive
<u>Buffer (R1a &amp; R2b)</u> Risk phrases:	R35; Causes severe burns
Safety phrases:	S24/25; Avoid contact with eyes and skin S26; In case of contact with eyes, rinse immediately with plenty of water and seek medical advice S36/37/39; Wear suitable protective clothing, gloves and eye / face protection S45; In case of accident or you feel unwell seek medical advice immediately (show label where possible).
Precautions:	See Section 3

16. Other information	
	The information herein is believed to be correct as of the date hereof but is provided without warranty of any kind. The recipient of our products is responsible for observing any National laws and guidelines applicable.