




Qnostics

SARS-CoV-2 (COVID-19) Control



RANDOX

In response to the current Coronavirus outbreak, Qnostics, a world leading manufacturer of QC solutions for molecular infectious disease testing has developed a SARS-CoV-2 control designed for use alongside molecular assays used in the testing of COVID-19. Accurate and reliable laboratory testing is an essential aspect of COVID-19 disease management and outbreak control. The use of a third-party quality control will allow diagnostic laboratories to effectively monitor the performance of diagnostic methods.

Why Choose Qnostics?

Whole Pathogen Control

Contains the entire SARS CoV-2 genome meaning it is compatible with the majority of commercial & in-house assays targeting the CDC and WHO consensus sequences.

Full Process Control

Designed to mimic the patient sample, full process controls will effectively monitor the performance of the entire testing process from extraction to amplification and detection.

Highly Characterised

Controls quantified by digital PCR ensuring batch to batch reproducibility.

True Third Party

An independent, unbiased assessment of assay performance is ensured in line with ISO 15189:2012 regulatory requirements.

Superior Manufacturing

Manufactured to ISO 13485:2016 standard and traceable to an internal reference preparation ensuring metrological traceability of test results.

Liquid for Ease-of-use

Conveniently supplied in a liquid frozen format meaning there is no additional preparation or handling required.

Negative Control Available

Complementary negative controls available offering a complete testing package.

SARS-CoV-2 Medium Q Control

Target Pathogen – SARS-CoV-2 (COVID-19)

Target – Whole Pathogen

Matrix – Transport Medium

Format – Liquid Frozen

Stability – Single use. Once thawed use immediately.

Regulatory Status – RUO

CATALOGUE NUMBER	PRODUCT DESCRIPTION	PACK SIZE
SCV2MQC	SARS-CoV-2 Medium Q Control	5 x 1 ml
TMNQC	Negative Q Control	5 x 1 ml

RANDOX

Copyright © 2019 Randox Laboratories Ltd. All rights Reserved. VAT number: GB 151682708.
Product availability may vary from country to country. Some products may be for Research Use Only.
For more information on product application and availability, please contact your local Randox Representative.

marketing@randox.com | randox.com