

|   |   |  |
|---|---|--|
| <b>Product</b>                              | : | AMOX – HRP Conjugate   |
| <b>Cat. No.</b>                             | : | HRP9355  |
| <b>Lot No.</b>                              | : | Lot dependent  |
| <b>Volume</b>                               | : | 0.5 mL   |
| <b>Recommended Working Dilution (ELISA)</b> | : | 1/160 when used in a competitive ELISA with polyclonal antibody PAS9995/PAS9996. |

The conjugate is supplied as a concentrate.  
Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.

|                |   |  |
|----------------|---|--|
| <b>Storage</b> | : | This conjugate may be stored for up to 3 months at <b>+2 - +8°C</b> .<br>For long term storage, aliquot and store at <b>≤-20°C</b> . |
|----------------|---|--|

*Product should be protected from light exposure.  
Avoid repeated freeze/thaw cycles.  
Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.*

**Notes & Precautions:**

This product as supplied is intended for research applications only, not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox BioReagents. A safety data sheet (SDS) can be supplied upon request.

26-May-22  
Information correct at time of going to print