

Product	:	Dinitrosalicylic Acid Hydrazide (DNSH) –
Cat. No.	:	HRP9396
Lot No.	:	Lot dependent
Volume	:	0.5 mL
Recommended Working Dilution (ELISA)	:	1/640 when used in a competitive ELISA with polyclonal antibody PAS9455/PAS9454.

The conjugate is supplied as a concentrate.
Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.

Storage	:	This conjugate may be stored for up to 3 months at +2 - +8°C . For long term storage, aliquot and store at ≤-20°C .
----------------	---	--

*Product should be protected from light exposure.
Avoid repeated freeze/thaw cycles.
Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.*

Notes & Precautions:

This product as supplied is intended for research applications only, not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox BioReagents. A safety data sheet (SDS) can be supplied upon request.

26-May-22
Information correct at time of going to print