

Interleukin-10

Human, Recombinant, *E. coli*

Cat. No.:	RCP9070	Size:	10µg
	RCP9071		100µg
	RCP9156		1000µg


Synonym: IL-10; IL10; Cytokine synthesis inhibitory factor; CSIF.

Description: Interleukin-10 is a multifunctional cytokine produced by a variety of cell types including T helper cells activated B cells and activated macrophages. Interleukin-10 regulates immune-mediated inflammation and modulates the function of cells such as lymphocytes, monocytes, natural killer cells and dendritic cells.

RANDOX recombinant Interleukin-10 comprises a 160 amino acid fragment (19-178) corresponding to the mature Interleukin-10 protein and is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

Form: Liquid

Purity: >95% by SDS-PAGE



References: Fickenscher, H. *et al.* (2002) *Trends Immunol.* **23(2)** : 89-96.