

---

**L-Selectin**

**Human, Recombinant, *E. coli***

**Cat. No.:** RCP9251

**Size:** 10µg

---

**Synonym:** Leukocyte Selectin, L-SEL, Leukocyte adhesion molecule LAM-1, Leukocyte cell adhesion molecule LECAM-1, TQ1, CD62L, Leukocyte surface antigen Leu-8, DREG, lymph node homing receptor LNHR, MEL-14

**Description:** L-Selectin belongs to a family of divalent cation-dependent carbohydrate-binding glycoproteins or adhesion molecules. L-Selectin is expressed constitutively on lymphocytes, monocytes and granulocytes and interacts specifically with carbohydrate groups on activated endothelial cells. L-Selectin may be shed by proteolytic cleavage and circulating levels in biological fluids may be used as an indicator of various pathological conditions.

**RANDOX** recombinant human L-selectin comprises a 294 amino acid fragment (39-332) corresponding to the mature L-Selectin extracellular domain and is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

**Form:** Liquid.  
Supplied in 1x Laemmli Buffer (25mM Tris-HCl pH6.8, 50mM DTT, 1% (w/v) SDS, 0.1% (w/v) Bromophenol Blue, 2.5% Glycerol).

**Purity:** Single band on Western blot.

**References:** Tedder, T.F *et al.* (1990) *J. Immunol.* **144** (2): 532-540.  
Schleiffenbaum, B. *et al.* (1992) *J. Cell Biol.* **119** (1): 229-238.