

MIP-1 α

Human, Recombinant, *E. coli*

Cat. No.:	RCP9022	Size:	10 μ g
	RCP9023		100 μ g
	RCP9140		1000 μ g


Synonym: Macrophage inflammatory protein 1-alpha; Chemokine (C-C motif) ligand 3; CCL-3; Small inducible cytokine-3; SCAY3; Tonsillar lymphocyte LD78 alpha protein; G0/G1 switch regulatory protein 19-1; SIS-beta; GOS19-1 protein; PAT464.1.

Description: MIP-1 α is a member of the CC chemokine family of proteins which have been shown to be produced by numerous cell types including activated T cells, B cells, monocytes, mast cells and endothelial cells. MIP-1 α acts as a monocyte chemoattractant and has been shown to act as a stem cell inhibitory factor.

RANDOX recombinant MIP-1 α comprises a 66 amino acid fragment (27-92) mature MIP-1 α protein and is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

Form: Liquid

Purity: >95% by SDS-PAGE



Reference(s): Menten, P. *et al.*, (2002) *Cytokine Growth Factor Rev.* 13 (6) : 455-481.