

Product : Tramadol – HRP9303 HRP Conjugate

Cat. No. : HRP9303

Lot No. : Lot dependent

Volume : 0.5 mL

Recommended Working Dilution (ELISA) : 1/320 when used in a competitive ELISA with polyclonal antibody PAS9293/PAS9278.

The conjugate is supplied as a concentrate.
Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.

Storage : This conjugate may be stored for up to 3 months at **+2 - +8°C**.
For long term storage, aliquot and store at **≤-20°C**.

*Product should be protected from light exposure.
Avoid repeated freeze/thaw cycles.
Avoid using azide containing buffers.
Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.*

Notes & Precautions:

This product as supplied is intended for research applications only, not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox BioReagents. A safety data sheet (SDS) can be supplied upon request.

30-May-22

Information correct at time of going to print