

| | | |
|---|---|---|
| Product | : | Triamcinolone Acetonide – HRP Conjugate |
| Cat. No. | : | HRP9501 |
| Lot No. | : | Lot dependent |
| Volume | : | 0.5 mL |
| Recommended Working Dilution (ELISA) | : | 1/2560 when used in a competitive ELISA with polyclonal antibody PAS10138/PAS10139. |

The conjugate is supplied as a concentrate.
Only dilute the concentrate as required and use the working strength conjugate immediately after dilution.

| | | |
|----------------|---|--|
| Storage | : | This conjugate may be stored for up to 3 months at +2 - +8°C . For long term storage, aliquot and store at ≤-20°C . |
|----------------|---|--|

*Product should be protected from light exposure.
Avoid repeated freeze/thaw cycles.
Avoid using azide containing buffers. Vial should be centrifuged briefly before opening to ensure all material is removed from the vial cap.*

Notes & Precautions:

This product as supplied is intended for research applications only, not for use in therapeutic or diagnostic applications without the expressed written authorization of Randox BioReagents. A safety data sheet (SDS) can be supplied upon request.

30-May-22
Information correct at time of going to print