

Eotaxin

Human, Recombinant, *E. coli*

Cat. No.:	RCP9049 RCP9050 RCP9149	Size:	10µg 100µg 1000µg
------------------	-------------------------------	--------------	-------------------------


Synonym: CCL11; Chemokine (C-C motif) ligand 11; Eosinophil chemotactic protein; Small inducible cytokine subfamily A member 11; SCYA11; Eotaxin precursor, Eosinophil chemotactic protein, Eot.

Description: Eotaxin / CCL11 is a member of the beta (C-C) motif chemokine family of inflammatory and immunoregulatory cytokines. Eotaxin is a potent and selective chemoattractant for eosinophils, but not monocytes, lymphocytes or neutrophils and is primarily produced by epithelial and phagocytic cells.

RANDOX recombinant Eotaxin comprises a 74 amino acid fragment (24-97) corresponding to the mature Eotaxin protein and is expressed in *E. coli* with an amino-terminal hexahistidine tag. This product is for research use only and is not intended for diagnostic or therapeutic use.

Form: Liquid

Purity: >95% by SDS-PAGE



References: Rankin, S.M. *et al.* (2000) *Mol. Med. Today* **6** (1): 20-27.